

## Immortalization of MEF with SV40 T antigen

---

### Reagents:

### Medium:

DMEM  
10% FCS  
1X Pen/strep  
1X Glutamine  
1X Non-essential amino acids  
55uM beta-mercaptoethanol

### Tag expression vector

Fugene

Serum free dMEM

---

### Method:

- Split primary MEF at pass 2-4 into 6 well dishes.
- Make 2 dilutions ~1/4 and 1/6 for each genotype.
- Grow overnight in 37C incubator
- Transfect ~25% confluent well with 2ug Tag expression vector using Fugene:
  - Change medium 1 hr before transfection. Use 2ml for each well.
  - Mix 2ug Tag expression vector with 20ul serum free dMEM in 1.5 ml tube.  
(Make a mix if you will be transfecting cells from several mice.)
  - In a separate 1.5ml tube, place 200ul serum free dMEM and then add 6ul Fugene (prewarmed to room temperature for 15 min).
  - Add DNA to Fugene mix and let stand for 15min at room temp.
  - Add mix dropwise to cells.
- Incubate cells 14-16hrs (overnight) in incubator.
- Change medium.
- When cells are just confluent (~2 days after transfection) split them into a 10cm dish=P1.

Additional splitting: The goal is to get cells that have been split 1/10 at least 5 times (i.e. a 1/100,000 fold splitting of the original cells) but if a low number of cells were transfected in the initial experiment then there may be few surviving cells from the early 1/10 splits. Because of this it is best to make a low density and high density plate at each pass. If it is clear that the low density (1/10) pass is growing, then the high density (1/4) plate is frozen as a backup or discarded. Often this 1/4 plate must be split before you can tell if the 1/10 plate is going to grow. (The P2-P3 low density plates can take 3-10 days to grow.)

\*\*\*\*\*Note that the 1/10 split is a strong negative selection against non-transformed cells whereas non-transformed cells can survive many 1/4 splits so that even if you keep the

high density plates as a back up, you will only get rid of the non-transformed cells by multiple rounds of 1/10 splitting. About 85% of the time I can use straight 1/10 splitting, but a few times I have used the 1/4 plate from pass 2 or 3 followed by 4-5 more passes of 1/10. By P6-100,000 (see below) the cells should grow to confluence in ~3-4 days after a 1/10 split.

